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Preparation of tooth surfaces after removal of fixed appliances

The adhesive removers are ideal for very carefully removing adhesive residue from the tooth surface. The non-cutting tip prevents injury to the gingival margin.

Both instruments are available in FG and RA shank.

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Titanium-nitride coated tungsten carbide cutters ensure very efficient preparation.

The titanium-nitride coating gives the instrument a smooth surface, which prevents the cutter clogging during preparation.

T23RX-012-FG		opt. 160.000 rpm
T21RX-012-FG		opt. 160.000 rpm
T21RX-012-RA		opt. 30.000 rpm
T21RX-010-RA		opt. 30.000 rpm

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The tungsten carbide finishers are for removing very small particles of adhesive residue.

H284-014-FG		opt. 160.000 rpm
H283K-016-FG		opt. 160.000 rpm
H282-010-RA		opt. 30.000 rpm
H282-012-RA		opt. 30.000 rpm

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The prepolisher (pink) from the DiaGloss range can be used as an alternative to tungsten carbide finishers. This instrument also effectively removes very small particles of adhesive residue without uncomfortable heat build-up.

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Polishing following removal of adhesive residue

The DiaGloss polisher produces a perfect high-lustre on the tooth surface. There is no uncomfortable heat build-up, as with conventional polishers.

The P19032, P19035 and P19039 polishers in particular are recommended from the large variety of shapes in the DiaGloss polisher range.

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CeraGlaze P30032

The very fine abrasive property of the polisher smooths scores and scratches left after removal of adhesive residue from the tooth surface.

The diamond grit produces a natural lustre on the enamel. The polishers can be used either wet or dry.

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
*Or-the-dens kit designed
by Dr Knut Thedens for removing
bracket adhesive.*

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

Rounding off • Stripping

The diamond strips with a saw are used in the case of slight over-crowding of the teeth, which prevents perfect rounding off of the dental arch and which has to be remedied by stripping.

The finer diamond coating of the FS3-FX ensures that less tooth structure is removed. The saw integrated in the diamond strips facilitates interdental separation.

In addition to the proven FS3-X there are now FS 6 diamond strips available in the same grit size, but without the integrated saw. The double width of the FS6 diamond strips ensures even more rapid separation.

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
*Perforated diamond finishing strip
with and without serrated edge.*

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
P6820

The brush is intended for cleaning rotary instruments. The brush also aids more effective cleaning of the instruments without damaging them, even with stubborn residue.

Shank	ISO	1
REF	P6820	

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Disinfection

Alcoholic-free solutions should be used for disinfection, as these do not damage the polishers and tungsten carbide instruments.

We highly recommend the Metasys cleaning system.

For more details for disinfection and sterilization refer to our catalogue „Rotary Dental Instruments Dental Surgery“.

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Wooden bur block

The high-quality, crafted, beech wood bur block is ideal for cleanly storing up to 25 instruments with HP shank.

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
© NTI-Kahla GmbH Rotary Dental Instruments • Printed in Germany

All rights reserved. No part of this catalogue may be reproduced in any form without written permission from the Management of NTI-Kahla GmbH. Reprint, translation, storing and data processing are not permitted without prior approval.

Products and colours may be subject to alterations. Printing errors excepted.

KFO - Labor • Orthodontic Laboratory


- L4 Gipsbearbeitung
Preparation of plaster/stone

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


- L11 Tiefgezogene Schienen
SoftWizard
*Vacuum-formed bite raisers
SoftWizard*

- L12 Soft Pol in drei abrasiven Stufen
SoftPol, available in three levels of abrasive

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- Diamond-coated dressing stone
Polishing of acrylic*

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- L15 *Preparation of metal
Polishing of metal*

KFO - Labor • Orthodontic Laboratory

Gipsbearbeitung • Preparation of plaster/stone

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Preparation of plaster/stone

The coarse cross-cut with additional cross blades ensures high removal with a very smooth surface when preparing Class 1 to 4 dental stones. The cross blades also prevent edge fracturing, which also makes the instrument ideal for fabricating study models.

A new more slimline cutter, HF079GE-045, with the same cross cut has been developed for also allowing preparation in areas that are unsuitable for large plaster cutters. This cutter can be used for perfect trimming of mucobuccal folds and labial frenums from stone models without fracturing the edges.

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Preparation of plaster/stone

A feature of these cutters is their amazing smooth running despite their 'martially' looking toothing.

The cutter does not vibrate and produces a very smooth surface.

Adhere to the recommended speed!

HF072SXCN-060  1  opt. 15.000 rpm

HF251SXCN-060  1  opt. 15.000 rpm

HF351SXCN-070  1  opt. 15.000 rpm

KFO - Labor • Orthodontic Laboratory

Gipsbearbeitung • Preparation of plaster/stone

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Set-Up

806.104.**365**.524.450

This diamond disc is impressive due to its perfect, smooth running and is ideal for basal separation of models.

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
806.104.**400**.514.220

A thickness of only 0.15 mm ensures minimal loss of structure when separating stone models interdentally.

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The large number of perforations provide an excellent view when preparing the surface of the model to be cut. It is also suitable for separation of Schwarz expansion plates.

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
806.104.**705**.514.220

This diamond disc has the same characteristics as the above mentioned diamond disc, but without any perforations.

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
806.104.**605**.514.220

The (non-toothed) perforated diamond disc is available as an alternative to toothed diamond discs. The fully diamond-coated disc has a thickness of 0.15 mm.

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Separating and trimming of wires and solder joints

Fibre-reinforced discs

The SD7000, which is only 0.2 mm thick, is ideal for separating custom-fabricated protrusive bars with bite jumping appliances. The SD7003 (0.5 mm) and SD7005 (1.0 mm) are used, e.g. for trimming solder joints or rounding off the tips of wires.

All three discs are fibre-reinforced and therefore guarantee a very high stability combined with minimal risk of fracture.

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Preparation of acrylic

The NTI Millennium Cutters in standard CE cross cut are recommended for preparing orthodontic acrylics cured in a pressure pot.

The extreme surface hardness of 3,600 HV ensures a long service life and is highly cost-effective. It also prevents edges fracturing during the preparation of acrylic.

A special production process guarantees that the cutters are very smooth-running, ensuring that the wrist is protected during preparation and preventing pain in the joints.

The shapes MC138CE-016, MC257CE-023, MC293CE-023, MC295CE-023 and MC261CE-023 are ideal for preparing areas difficult to access.

The FE cut is used for efficient preparation of acrylic.

Advantage:

The FE cut removes less material than when using the CE cut and also produces a very fine surface.

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HF257CE-060/HF257CNR-060

Both cutters are ideal for gentle preparation of the junction between acrylic and wire elements, depending on the type of surface required. The cutters operate less aggressively in the upper quarter of the working section. This considerably reduces the risk of inadvertently damaging wire elements and consequently producing break-off points.

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Preparation of acrylic

Special blades for gentle, smooth adjusting of denture acrylics.

The newly developed blades ensure quick, economic removal of material.

The blade design prevents the usual white excess with acrylic adjusting.

This cutter is very quiet and removes material efficiently without leaving flash. A smooth surface remains as a result, which simplifies subsequent polishing.

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U - BOW

The repair cutter instrument is suitable for placing a U-bow in the U-bow activator developed by Prof. Karwetzky.

This cutter facilitates preparing hollows in areas of the acrylic where the U-bow is to be incorporated at a later stage.

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Bite surfaces

Acrylic Grinder
The diamond-coated hollow cutter is recommended for contouring an anterior acrylic bite-raiser.

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
806.104.353.524.220

806.104.396.514.220

The fully diamond-coated discs are perforated, which facilitates precision preparation.

806.104.362.524.100

The diamond disc 806.104.362.524.100 is used for areas difficult to prepare, as the special diamond coating of the disc ensures cooler preparation and prevents “sticking” due to heat build-up.

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Permanently soft acrylics/ soft relines

Grinding caps are the most efficient instruments for preparing positioners fabricated from materials applied from a cartridge. Grinding caps are not only excellent for material removal, but they also produce a uniform surface.

The use of cutters with the QSCN coarse straight blade cross cut produces very good results if the cutters are initially heated to “operating temperature”. This can be achieved by trimming standard hard acrylic with the cutter for approx. one minute at a speed of 15000 – 20000 rpm.

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
Vacuum-formed bite raisers

“Hard bite raisers” are easily prepared with the Fissure Bur.

Edges on the margins of the bite raiser do not necessarily require further preparation, as the Fissure Bur produces a very smooth, uniform surface.

The HF515S-023 instrument is also used with bite raisers fabricated using a combination of hard and soft acrylics, e.g. Durasoft by Scheu-Dental GmbH.

When trimming splint peripheries, flash is often produced on the soft acrylic, which up until now has been difficult to remove. Using the G9006 the peripheries of the splint, including the interdental areas, can be prepared to produce a very smooth surface. The SuperMax G9006 is a sinter diamond trimmer with organic binders, which absorb the heat during trimming. The speed should never exceed 10000 rpm. It is important always to trim from the soft acrylic towards the hard acrylic. Low pressure application combined with a low speed ensures a high service life of the instrument.

SoftWizard

Flash created when preparing soft acrylics can be easily removed using the SoftWizard P2100.

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
P4203

Vacuum-formed splints

This cutter is ideal for orthodontic acrylic and vacuum-formed splints; it produces a polished-like surface.

The cutter is very quiet running.

It is also ideal for smoothing the ends of wire. Pressure points can be removed at the chairside with out having to re-polish.

The HF138NP-023, in particular, is ideal for accessing very fine niches.

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SoftPol

G: ideal for pre-contouring, effectively and quietly.

M: for pre-smoothing.

F: preparation without feather edges or streaks on the surface, produces a high-quality, smooth surface.

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Diamond-coated dressing stone

If the SoftWizard is tapered using the P4060 diamond-coated dressing stone, it can also be used for efficiently preparing interdental areas and producing a very smooth surface. Grinding caps are used for preparing bite raisers fabricated using permanently soft acrylic.

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Polishing of acrylic

NTI AcrylicMaster

The NTI AcrylicMaster P0634 and P0644 (two-stage polishing system) and NTI Prosthetics P0674, P0664 and P0654 (three-stage polishing system) produce a high removal performance in minimum time when polishing acrylic. The acrylic is prepared in sequence from coarse to fine.

In particular the P0674, P0664, P0654, P0679, P0669 and P0659 polishers from the large variety of shapes in the NTI Prosthetics range ensure highly efficient preparation due to their slimline design.

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
Preparation of metal

The SuperMax is ideal for rounding off the cutting edges on wires. Particular features of the SuperMax are its excellent material removal with smooth ground surfaces and its outstanding service life.

The G9006, G9007 and G9008 instruments from the large variety of shapes in the SuperMax range are ideal for preparing large surface areas.

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The superfine cross-cut Millennium cutter can be used for initial preparation of laboratory fabricated metal appliances (Quad-helix /QH and Crozat appliances), which still require soldering.

The new shape of the instrument ensures easy preparation, even in areas difficult to access. Outstanding material removal with a smooth surface facilitates subsequent polishing.

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Polishing of metal

The two-stage NTI CeraGlaze is ideal for producing a perfect high-lustre polish at the solder joints. Another feature of the NTI CeraGlaze is its extremely long service life.


The prepolisher P3044, which is slightly abrasive, produces a smooth, matt sheen in minimum time.

The P30044 produces a mirror smooth surface. Further polishing is not required.

The P3042, P30042, P3041, P30041, P3001 and P30001 polishers in the comprehensive NTI CeraGlaze range also ensure perfect surface preparation.

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
This brush is an ideal aid in the laboratory. Cleaning cutters with this brush increases their service life.

L15

Shank  
 **P6820**