


Not Only a **new Class II** **Device**



Features


- Simple
- Precise
- Adaptable
- Efficient



Introduction
by Dr.Yoav Mazor

MULTI-FA

Multi functional Appliance
for CL II and CL III
Correction



Experts in Orthodontics

GC Orthodontics Europe GmbH

www.gc.dental/ortho

MULTI-FA

Double hinge for an anatomical positioning without compromise



Improvement of the biomechanical control



Fix the hinge upon treatment plan requirements



One appliance for 90% cases




mesio – distal stop made with a flow



or a crimpable stop

MULTI-FA for tubes, adaptable on patients



92-5001-0000 MULTI FA REGULAR (left & right)

92-5002-0000 MULTI FA SHORT (left & right)

92-5003-0000 MULTI FA LONG (left & right)

92-5004-0000 MULTI FA BAND-R (left & right)

92-5005-0000 MULTI FA CRIMPABLE STOP (x10)

